

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com

Distal Clavicular Locking Plate

<i>Products</i>	<i>Picture</i>	<i>REF</i>	<i>Specification</i>	<i>Thickness</i>	<i>Width</i>	<i>Length</i>
Distal Clavicular Locking Plate (Use 2.7/3. 5 Locking Screw/3.5 Cortical Screw)		01.1158-01	3 holes L	4.2	11	61
		01.1158-02	4 holes L	4.2	11	71
		01.1158-03	5 holes L	4.2	11	81
		01.1158-11	3 holes R	4.2	11	61
		01.1158-12	4 holes R	4.2	11	71
		01.1158-13	5 holes R	4.2	11	81

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com
