

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com

Calcaneous Locking Plate

<i>Products</i>	<i>Picture</i>	<i>REF</i>	<i>Specification</i>	<i>Thickness</i>	<i>Width</i>	<i>Length</i>
Calcaneous Locking Plate(Use 3.5 Locking Screw)		5100-4905	Short 3 holes Left	1.5	8.5	76.5
		5100-4906	Short 3 holes Right	1.5	8.5	76.5
		5100-4907	Long 3 holes Left	1.5	8.5	76.5
		5100-4908	Long 3 holes Right	1.5	8.5	76.5
		5100-4909	Short 4 holes Left	1.5	8.5	86.5
		5100-4910	Short 4 holes Right	1.5	8.5	86.5
		5100-4911	Long 4 holes Left	1.5	8.5	86.5
		5100-4912	Long 4 holes Right	1.5	8.5	86.5

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com
