

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com

S-clavicle Locking Plate-III

<i>Products</i>	<i>Picture</i>	<i>REF</i>	<i>Specification</i>	<i>Thickness</i>	<i>Width</i>	<i>Length</i>
S-clavicle Locking Plate-III (Use 3.5 Locking Screw/3.5 Cortical Screw)		5100-8701	6 holes L	3.2	10	92
		5100-8702	7 holes L	3.2	10	105
		5100-8703	8 holes L	3.2	10	119
		5100-8704	6 holes R	3.2	10	92
		5100-8705	7 holes R	3.2	10	105
		5100-8706	8 holes R	3.2	10	119

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com
