

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com

Humeral External Fixator

<i>Products</i>	<i>Picture</i>	<i>REF</i>	<i>Specification</i>
Humeral External Fixator		9100-2101	Rod-rod Coupling 8/8mm 2PCS Pin-rod Coupling 8/4-5mm 6PCS Connecting Rod 8*250mm 2PCS Connecting Rod 8*200mm 1PCS Cortical Pin 5*130mm 6PCS

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com
