

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com

Knee Joint External Fixator

<i>Products</i>	<i>Picture</i>	<i>REF</i>	<i>Specification</i>
<i>Knee Joint External Fixator</i>		<i>9100-2701</i>	<i>Rod-rod Coupling 8/8mm 8PCS</i> <i>Pin-rod Coupling 8/4-5mm 14PCS</i> <i>Semi-circular Curved Rod 8mm Large 2PCS</i> <i>Curved Rod 8mm 2PC</i> <i>Spanning Joint Connecting Rod 8mm 2PCS</i> <i>Cortical Pin 5*130mm 14PCS</i>

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com